

RIBCARE

TIR14A1



RibCare

The RibCare System consists of locking plates, locking screws for the fixation and stabilization of ribs. The long wand and flexible driver is developed for the surgical treatment of posterior rib specifically. The flexible driver is patented as a surgical instrument(10-2017-0115351).

Indications


The RibCare System is indicated for the fixation and stabilization of rib fractures, fusions and osteotomies of bone.

Features

- All implants are manufactured from titanium alloy(grade 23) for flexibility, strength, biological feasibility.
- The long wand and the Flexible driver can reach to posterior part of the rib.
- The screw has no need pre-drill process for fixation of the plate.

RibCare System

RibCare Plate



Part No	길이(L)mm	홀수N
TIR14A1-PP1	38	5
TIR14A1-PP2	53	7
TIR14A1-PP3	68	9
TIR14A1-PP4	98	13
TIR14A1-PP5	128	17

※ TIR14A1-PP1~PP5의 형상은 모두 동일하며 제품의 전체 길이와 홀의 개수만 다름.


RibCare Screw



Part No	길이(L)mm	직경(Ø)mm
TIR14A1-PS0	6	3.4
TIR14A1-PS1	8	3.4
TIR14A1-PS2	10	3.4
TIR14A1-PS3	12	3.4
TIR14A1-PS4	14	3.4
TIR14A1-PS5	16	3.4
TIR14A1-PS6	18	3.4
TIR14A1-PS7	20	3.4


※ TIR14A1-PS1~PS7의 형상은 모두 동일하며 제품의 전체 길이만 다름.

RibCare Sternal Plate



Part No	길이(L)mm
TIR14A1-SL1	62
TIR14A1-SL2	81
TIR14A1-SL3	100
TIR14A1-SS1	28.7

RibCare Mini Tray



Component

RibcarePlate	RibCareScrew
TIR14A1-PS0	TIR14A1-PP1
TIR14A1-PS1	TIR14A1-PP2
TIR14A1-PS2	TIR14A1-PP3
TIR14A1-PS3	TIR14A1-PP4
TIR14A1-PS4	TIR14A1-PP5
TIR14A1-PS5	

RibCare Instrument



RibCare Tray



TIR14A1-FD4 -
MINIPLATE WRENCH DRIVER



TIR14A1-ITA -
T DRIVER




TIR14A1-FDA - FLEXIBLE DRIVER A TYPE



TIR14A1-FDB - FLEXIBLE DRIVER B TYPE



TIR14A1-ISA - MINIPLATE SHORT WAND



TIR14A1-IRG - RIBGAGE



TIR14A1-ILA - MINIPLATE LONG WAND



TIR14A2 - POWER DRIVER



TIR14A2-PD6 - BOX DRIVER




TIR14A2-PD7 - CROSS DRIVER



TIR14A1-IBA - MINIPLATE BENDER

RIBCARE_{TIR14A1}



Distributor

 Hong-Eun Medical Co.,Ltd

Tel : +82-2-455-4011

www.hongunmedical.co.kr

5F, Hong-eun Bldg, 176, Jayang-ro

Gwangjin-gu, Seoul,05038, Korea